



Cleanroom Swabs

Microdenier | Polyester | Foam | ESD | Cotton | Specialty

Sterile | Sterile Dry Collection and Transport System | Cleaning Validation & TOC Validation Kits Alpha[®] | Absorbond[®] | CleanFoam[®] | ESD-Safe | Cotton | Specialty | General Purpose

Table of Contents

Sterile
Sterile Dry Collection and Transport System Pg. 3
Cleaning Validation Pg. 4
Cleaning Validation Kits Pg. 4
Microdenier
Polyester
CleanFoam [®]
ESD-Safe
Cotton
Specialty Pg. 10
General Purpose Pg. 11

Our highly professional team will help you select and develop products for any critical environment application. Contact information is located on the back of this brochure or visit us at **www.texwipe.com**.

About Texwipe

Texwipe's commitment to innovation, leadership and quality in cleanroom consumable products spans more than fifty years. We invest in technology to respond to our customers' evolving needs in contamination control.

Texwipe's Core Values are:

- **Innovation** Texwipe pioneers the latest technologies to provide innovation in contamination control products and processes.
- Quality Texwipe product quality is maintained by the most advanced testing and quality control standards in the industry.
- Technology Leadership Texwipe leads our industry in testing metrology, methods and processes to reduce contamination.

Throughout Texwipe's global operations, we support our customers with products designed to exceed the requirements for cleanroom consumable products. Our highly professional team will help you select and develop products for any critical environment application.

Texwipe CleanTips®

Texwipe Cleanroom Swabs, with their trademark green handles, are the standard in cleanrooms across the globe. Engineered for superior cleanliness and consistency for a variety of applications, Texwipe has a swab that will provide the best contamination performance available.

Our products range from Microdenier, Polyester, CleanFoam[®], ESD-Safe, General Purpose, Cotton and our Specialty line featuring the pre-wetted CrushTube[™] Swab.

Texwipe Cleanroom Swabs are manufactured with a complete thermal bond, eliminating adhesive contamination*. Our swabs are manufactured using high-precision automated processes, resulting in swabs constructed to exacting and consistent tolerances, lot coded for traceability and quality control.

Custom Swabs

At Texwipe[®] we understand that the best swab for your application may be something unique. We can create a custom swab for your specific applications.

Sterile Swabs

We can create sterile versions of all Texwipe swabs or something unique and specific to your application.

*excludes Cotton Swabs



www.texwipe.com



Sterile Swabs

Texwipe's Sterile Swabs are available in USP-grade cotton, polyurethane foam and medical-grade spun polyester. The swabs are made to exacting and consistent tolerances using high-precision automated processes. The easy-to-open, individually wrapped swabs ensure sterility at point-of-use. The swabs are visible through clear packaging which also offers visibility of seal integrity.



 $\langle \neg \rangle$

Cotton Series

Part No. STX705P	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	7.0mm (0.276")	7.0mm (0.276'')	18.2mm (0.717")	2.5mm (0.098")	2.5mm (0.098'')	134.0 mm (5.276'')	155.0mm (6.102'')	1 swab per sleeve 50 swabs per box 10 boxes per case
Part No. STX705W	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging

Foam Series

Part No. STX708A	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	19.7mm	9.4mm	26.4mm	6.6mm	3.2mm	105.8mm	132.2mm	1 swab per sleeve
	(0.776")	(0.370")	(1.040")	(0.260'')	(0.126")	(4.165")	(5.205")	50 swabs per box
								10 boxes per case
07/7408				I				
Part No. STX712A	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm	7.8mm	25.7mm	5.2mm	3.0mm	101.8mm	127.5mm	1 swab per sleeve
	(0.500'')	(0.307")	(1.012")	(0.205")	(0.118")	(4.008")	(5.020'')	50 swabs per box
								10 boxes per case
							EXVIRE:	

Polyester Series

	5.8mm (0.228'')	5.8mm (0.228'')	18.0mm (0.709'')	2.5mm (0.098'')	2.5mm (0.098'')	134.0mm (5.276'')	155.0mm (6.102'')	1 swab per sleeve 50 swabs per box
	(0.220)	(0.220)	(0.705)	(0.050)	(0.030)	(5.270)	(0.102)	10 boxes per case
Part No. STX764	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
Part No. STX764	Head Width 4.6mm	Head Thickness 4.6mm	Head Length 17.0mm	Handle Width 2.5mm	Handle Thickness 2.5mm	Handle Length	Total Swab Length 155.0mm	Packaging 1 swab per sleeve
Part No. STX764			Ŭ			Ŭ	, in the second s	

Sterile Dry Collection and Transport System

Texwipe's Sterile Dry Collection and Transport System Swabs consist of sterile ready-to-use units for collection and transport of bacteriological specimens. The system is designed to maintain integrity of both sample collection and sample transport to the lab. Swabs are constructed from 100% USP-grade spun cotton or 100% USP-grade spun polyester. The swabs are bonded securely to a polystyrene shaft by an aqueous-based adhesive. 6.1" (155 mm) long swabs enable acess to hard-to-reach sampling sites. The swabs are made to exacting and consistent tolerances using high-precision automated processes.



| Page 3

Cotton Series

Part No. STX705PT	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	7.0mm (0.276'')	7.0mm (0.276'')	18.2mm (0.717")	2.5mm (0.098'')	2.5mm (0.098'')	134.0 mm (5.276')	155.0mm (6.102'')	1 swab per tube 100 tubes per box 5 boxes per case
03/2	xwipe	Kernersville, N www.texwipe.		Americas +1 80 Europe +31 88 Asia +63 49 54	1307 410			

Polyester Series

Part No. STX764T	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	4.6mm (0.181'')	4.6mm (0.181")	17.0mm (0.669'')	2.5mm (0.098'')	2.5mm (0.098'')	134.0 mm (5.276'')	155.0mm (6.102'')	1 swab per tube 100 swabs per box 5 boxes per case
103/2018	xwipe	www.texwipt	NC 27284 USA 5.com	Americas +1.8/ Europe +31.88 Asia +63.49.54	1307 410			
1								
					- ANDRE			
					- take	705PT		
			ARRENT			A	and a second	The second se
							76	at 1
							A CONTRACTOR	
	Sterile Av	reileble				and the second sec	and a	
		in sterile upon	request.		Tex	xwipe is certifi	ed to ISO 9001 a	nd ISO 13485.

Cleaning Validation Swabs

Engineered to deliver consistent and accurate sample recovery for cleaning validation for both specific (HPLC, UV-Vis, IMS) and nonspecific (TOC) analytical methods. The swabs are effective with a wide range of diluents.

Alpha[®]

Alpha series swab additionally processed to $<50 \ \mu g/L \ TOC$ per swab.

Part No. TX714K	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500'')	4.2mm (0.165'')	25.7mm (1.012'')	5.2mm (0.205'')	3.0mm (0.118'')	101.8mm (4.008'')	127.5mm (5.020'')	20 swabs per bag 50 bags per case
						III	wos -)	
Part No. TX715	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500'')	4.2mm (0.165")	25.7mm (1.012")	5.2mm (0.205'')	3.0mm (0.118")	101.8mm (4.008")	127.5mm (5.020'')	100 swabs per bag 2 inner bags of 50 swabs 10 bags per case
Part No. TX761K	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.8mm (0.268'')	2.8mm (0.110'')	16.8mm (0.661'')	3.2mm (0.126'')	3.2mm (0.126'')	145.5mm (5.728'')	162.3mm (6.390'')	20 swabs per bag 50 bags per case
and the second s	_							19805

Absorbond[®]

Surfactant-free Absorbond Series swab for HPLC/IMS analysis.

Part No. TX716	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500'')	4.2mm (0.165")	25.7mm (1.012'')	5.2mm (0.205'')	3.0mm (0.118'')	101.8mm (4.008'')	127.5mm (5.020'')	20 swabs per bag 50 bags per case
Farmer and the	(0.500)	(0.105)	(1.012)	(0.205)	(0.110)	(4.008)	(5.020)	SU Days per case
The second second	_	_	_	_		56303		

TOC Cleaning Validation Kits

Texwipe's TOC Cleaning Validation Kits are designed to simplify sampling as part of a cleaning validation protocol. The swabs and vials are further cleaned using a proprietary processes to make these kits excellent for use in Total Organic Carbon (TOC) analysis.

The kits are an organized, convenient and secure method for transporting swabs and vials from storage area to sampling sites to laboratory with minimal chance of contamination. Packaged in cleanroom-compatible, recyclable polypropylene boxes.

Available with Texwipe Cleanroom Swab TX714K (TX3340 and TX3342) or TX761K (TX3343).

Product Number	Description	Packaging
TX3340	TOC Cleaning Validation Kit: 12 vials + 24 swabs (TX714K) + 12 labels	18 kits/case
TX3342	TOC Cleaning Validation Kit, Bulk: 72 vials + 144 swabs (TX714K) + 72 labels	1 kit/case
TX3343	TOC Cleaning Validation Kit: 12 vials + 24 swabs (TX761K) + 12 labels	18 kits/case

All Texwipe swabs available in sterile upon request.

TX3340

TOC Cleaning

Validation Kit

Polyester Swabs

Microdenier Series

Texwipe's Microdenier Swab Series is made from the cleanest microdenier material, 100% polyester, nylon-free. The complete thermal bond construction eliminates adhesive contamination.

Cleanroom laundered, providing ultra low levels of particles, NVRs (non volatile residues) and ions, and engineered to provide enhanced absorbency and quick wicking performance. Streak-free cleaning with scratch-sensitive surfaces and compatible with QAC (quaternary ammonium chloride) cleaners and disinfectants. 100% virgin polypropylene handle.

Cleanroom manufactured, made to exacting and consistent tolerances using high-precision automated processes. Lot coded for traceability and quality control.

Packaged in a reclosable, silicone-free bag.







Alpha® Series

Texwipe's Alpha[®] Swab Series is made from the cleanest 100% polyester knit materials with a double layer of fabric for enhanced absorbency. A complete thermal bond construction eliminates adhesive contamination. Cleanroom laundered and manufactured, providing ultra low levels of contaminants. Engineered to resist snagging and abrading. 100% virgin polypropylene handle. Made to exacting, consistent tolerances using high-precision automated processes. Packaged in a silicone-free, amide-free bag.

Part No. TX714A	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm (0.500")	4.2mm (0.165")	25.7mm (1.012'')	5.2mm (0.205'')	3.0mm (0.118'')	101.8mm (4.008'')	127.5mm (5.020'')	100 swabs per bag 2 inner bags of 50 swabs
							TET GEFALL	10 bags per case
Part No. TX743B	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	4.0mm (0.157")	3.6mm (0.142")	12.0mm (0.472")	3.0mm (0.118")	3.0mm (0.118'')	58.3mm (2.295'')	70.3mm (2.768'')	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. TX754B	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
_	3.5mm (0.138'')	3.0mm (0.118")	9.5mm (0.374'')	2.4mm (0.095")	2.4mm (0.095")	60.0mm (2.362'')	69.5mm (2.736")	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. TX758B	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
-	3.2mm (0.126'')	2.5mm (0.098'')	11.0mm (0.433")	2.2mm (0.087'')	2.2mm (0.087'')	60.0mm (2.362'')	70.0mm (2.756'')	500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. TX761	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.8mm (0.268'')	2.8mm (0.110'')	16.8mm (0.661'')	3.2mm (0.126'')	3.2mm (0.126")	145.5mm (5.728'')	162.3mm (6.390'')	100 swabs per bag 10 bags per case
			_				_	- SCHERC

Absorbond® Series

Texwipe's Absorbond[®] Swabs Series is made from 100% polyester (hydroentangled) nonwoven material double layer for enhanced absorbency. The complete thermal bond construction eliminates adhesive contamination. Cleanroom processed and manufactured, providing ultra low levels of contamination control. Engineered for use on sensitive surfaces. 100% virgin polypropylene handle. Made to exacting, consistent tolerances using high-precision automated processes. Packaged in a silicone-free, amide-free bag.

Part No. TX716	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	12.7mm	4.2mm	25.7mm	5.2mm	3.0mm	101.8mm	127.5mm	20 swabs per bag
and the Wild	(0.500'')	(0.165'')	(1.012'')	(0.205'')	(0.118'')	(4.008'')	(5.020'')	50 bags per case
Direct .	-			_	_	178	X///IRE	
State of the state								
Part No. TX759B	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.2mm	1.5mm	10.0mm	2.2mm	2.2mm	60.0mm	70.0mm	500 swabs per bag
	(0.126'')	(0.059'')	(0.394")	(0.087'')	(0.087")	(2.362")	(2.756")	5 inner bags of 100 swabs
								5 bags per case
				_				
Part No. TX762	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	6.8mm	2.6mm	16.8mm	3.2mm	3.2mm	145.5mm	162.3mm	100 swabs per bag
	(0.268'')	(0.102")	(0.661")	(0.126'')	(0.126'')	(5.728'')	(6.390'')	10 bags per case

CleanFoam[®] Swabs


CleanFoam® Series A

Texwipe's CleanFoam[®] Swab Series is made with 100% polyurethane foam, 100 pores per inch, for good chemical compatibility with a variety of solutions. The complete thermal bond construction eliminates adhesive contamination.

100% virgin polypropylene handle. Cleanroom processed and made to exacting, consistent tolerances using high-precision automated processes. Packaged in a silicone-free, amide-free bag.



Foam Series A features 100 swabs per bag.



CleanFoam® Swabs

CleanFoam® Series B

Texwipe's CleanFoam[®] Swab Series is made with 100% polyurethane foam, 100 pores per inch, for good chemical compatibility with a variety of solutions. The complete thermal bond construction eliminates adhesive contamination.

100% virgin polypropylene handle. Cleanroom processed and made to exacting, consistent tolerances using high-precision automated processes. Packaged in a silicone-free, amide-free bag.



۲.

Foam Series B features 500 swabs per bag.

art No. TX740B	Head Width 6.2mm (0.244")	Head Thickness 5.0mm (0.197")	Head Length 17.0mm (0.669")	Handle Width 3.2mm (0.126")	Handle Thickness 3.2mm (0.126")	Handle Length 146.0mm (5.748")	Total Swab Length 163.0mm (6.417")	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
art No. TX741B	Head Width 3.4mm (0.134")	Head Thickness 3.0mm (0.118")	Head Length 10.5mm (0.413")	Handle Width 2.4mm (0.094")	Handle Thickness 2.4mm (0.094")	Handle Length 56.8mm (2.236")	Total Swab Length 67.3mm (2.650")	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
irt No. TX742B	Head Width 3.6mm (0.142")	Head Thickness 3.5mm (0.138")	Head Length 12.0mm (0.472")	Handle Width 3.0mm (0.118")	Handle Thickness 3.0mm (0.118")	Handle Length 59.0mm (2.323")	Total Swab Length 71.0mm (2.795")	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
urt No. TX751B	Head Width 3.5mm (0.138")	Head Thickness 3.5mm (0.138")	Head Length 16.0mm (0.630")	Handle Width 2.4mm (0.094")	Handle Thickness 2.4mm (0.094")	Handle Length 51.0mm (2.008")	Total Swab Length 67.0mm (2.638")	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
rt No. TX752B	Head Width 3.4mm (0.134")	Head Thickness 3.4mm (0.134")	Head Length 20.0mm (0.787")	Handle Width 2.5mm (0.098")	Handle Thickness 2.5mm (0.098")	Handle Length 50.0mm (1.969")	Total Swab Length 70.0mm (2.756")	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
rt No. TX757B	Head Width 3.2mm (0.126")	Head Thickness 2.6mm (0.102")	Head Length 10.0mm (0.395")	Handle Width 2.2mm (0.087")	Handle Thickness 2.2mm (0.087")	Handle Length 59.7mm (2.350")	Total Swab Length 69.7mm (2.744")	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case

ESD-Safe Swabs

Texwipe's ESD Swab Series features a Stat-Rite[®] handle, a patented Inherently Dissipative Polymer, to provide clean ESD protection without contributing particulate contaminants to the process or environment. The complete thermal bond construction eliminates adhesive contamination. Packaged in an anti-static bag.

Available in a choice of CleanFoam®, Alpha® Polyester, or Stick.



CleanFoam® ESD

Part No. TX750E	Head Width 3.5mm (0.138'')	Head Thickness 3.5mm (0.138")	Head Length 10.0mm (0.394'')	Handle Width 3.0mm (0.118")	Handle Thickness 3.0mm (0.118")	Handle Length 78.5mm (3.090")	Total Swab Length 88.5mm (3.484")	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. TX753E	Head Width 3.7mm (0.146'')	Head Thickness 3.4mm (0.134")	Head Length 10.0mm (0.394")	Handle Width 3.0mm (0.118")	Handle Thickness 3.0mm (0.118")	Handle Length 80.0mm (3.150")	Total Swab Length 90.0mm (3.543")	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case
Part No. TX757E	Head Width 3.8mm (0.150'')	Head Thickness 3.4mm (0.134")	Head Length 10.0mm (0.394")	Handle Width 3.0mm (0.118'')	Handle Thickness 3.0mm (0.118")	Handle Length 80.7mm (3.177'')	Total Swab Length 90.7mm (3.571 '')	Packaging 500 swabs per bag 5 inner bags of 100 swabs 5 bags per case

Alpha® Polyester ESD

Part No. TX758E	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	4.0mm (0157'')	3.0mm (0.118'')	10.0mm (0.394'')	3.0mm (0.118'')	3.0mm (0.118'')	80.5mm (3.169'')	90.5mm (3.563'')	500 swabs per bag 5 inner bags of 100 swabs
								5 bags per case
Part No. TX761D	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
Part No. TX761D	Head Width 6.8mm	Head Thickness 4.0mm	Head Length 17.0mm	Handle Width 3.0mm	Handle Thickness 3.0mm	Handle Length	Total Swab Length 164.0mm	Packaging 100 swabs per bag
Part No. TX761D						Ŭ.	, in the second se	

Stick ESD

Part No. TX769E	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	1.0mm (0.039'')	1.0mm (0.039")	9.6mm (0.378'')	2.4mm (0.094'')	2.4mm (0.094'')	56.5mm (2.224'')	66.1mm (2.602'')	500 swabs per bag 5 inner bags of 100 swabs
								5 bags per case



Cotton Swabs

Texwipe's Cotton and Cotton/Foam Swab Series is designed to maximize the characteristics of cotton. Available in 100% cotton or cotton/polyurethane foam. These swabs are bonded securely to a wooden shaft by an aqueous-based adhesive.

 \Box

Part No. TX700B	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging	
	13.5mm (0.531'')	13.5mm (0.531'')	28.0mm (1.102'')	2.5mm (0.098'')	2.5mm (0.098'')	126.0mm (4.961'')	154.0mm (6.063'')	500 swabs per bag 10 inner bags of 50 swabs	
								5 bags per case	
	J.								
TV705	11	line of The Selence of	line of Lengeth	Handle Width	Usedle Thistory	llegalis Legentin	Table Quark Law atta	Destacion	
Part No. TX705	Head Width	Head Thickness	Head Length		Handle Thickness	Handle Length	Total Swab Length	Packaging	
	7.0mm (0.276'')	7.0mm (0.276'')	17.0mm (0.669'')	2.5mm (0.098'')	2.5mm (0.098'')	135.0mm (5.315'')	152.0mm (5.984'')	1,000 swabs per bag 10 inner bags of 100 swabs	
	(0.270)	(0.270)	(0.009)	(0.090)	(0.096)	(5.515)	(5.964)	5 bags per case	
Part No. TX705P	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging	
	7.0mm	7.0mm	17.0mm	2.5mm	2.5mm	135.0mm	152.0mm	1,000 swabs per bag	
	(0.276")	(0.276'')	(0.669'')	(0.098'')	(0.098'')	(5.315'')	(5.984'')	10 inner bags of 100 swabs 5 bags per case	
Part No. TX720B	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging	
	11.8mm	10.0mm	26.3mm	2.5mm	2.5mm	127.0mm	154.0mm	500 swabs per bag	
	(0.465")	(0.394'')	(1.035")	(0.098'')	(0.098'')	(5.000'')	(6.063'')	10 inner bags of 50 swabs 5 bags per case	

Specialty Swabs & Tools

TX726 is our innovative CrushTube Swab, featuring 91% IPA in an enclosed tube. Packaged in a resealable bag.

Part No. TX726	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	8.3mm	n/a	10.0mm	9.5mm	9.5mm	44.0mm	47.6mm	50 swabs per box
	(0.327")		(0.394")	(0.37")	(0.37")	(1.73")	(1.89")	10 boxes per case
. (.								

TX730 is a 100% nylon swab with micro bristles that grab particles or remove liquids by capillary forces. Ideal for tight spaces. Packaged in a silicone-free, amide-free bag.

Part No. TX730	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	0.8mm	0.8mm	15.5mm	2.8mm	2.8mm	152.5mm	168.0mm	50 swabs per bag
	(0.031'')	(0.031'')	(0.610'')	(0.110'')	(0.110'')	(6.004'')	(6.614'')	10 bags per case

Sterile Available All Texwipe swabs available in sterile upon request.

General Purpose Swabs

Texwipe's General Purpose Swab Series is designed for those customers needing a swab in industrial or general purpose applications. 100% polypropylene handle provides excellent chemical compatibility. Complete thermal bond construction eliminates adhesive contamination.

Packaged in a reclosable, silicone-free bag. Available in Polyester Honeycomb and Foam.

Polyester Honeycomb

Texwipe's General Purpose Polyester Swab Series, made with 100% polyester, was engineered with a honeycomb-type cell designed to maximize absorbency and capture of large particles and dirt in less critical areas.



Reclosable, silicone-free bag

Part No. TX801	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	13.5mm	3.5mm	25.2mm	5.2mm	3.0mm	102.3mm	127.5mm	100 swabs per bag
	(0.531")	(0.138'')	(0.992'')	(0.205'')	(0.118'')	(4.028")	(5.020")	6 bags per case
	-						_	
				_			-	
	-							

					Handle Length	Total Swab Length	i ackagilig
3.5mm	3.0mm	9.2mm	2.5mm	2.5mm	60.0mm	69.2mm	500 swabs per bag
(0.138'')	(0.118'')	(0.362")	(0.098")	(0.098'')	(2.362")	(2.724'')	5 bags per case



General Purpose Swabs



 \Box

Foam

Texwipe's General Purpose Foam Swab Series, made with 100% polyurethane foam, was engineered for applying and absorbing adhesives, greases and liquids in less critical areas.



Part No. TX804	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	3.5mm	3.5mm	15.5mm	2.4mm	2.4mm	51.0mm	66.5mm	500 swabs per bag
	(0.138'')	(0.138'')	(0.610'')	(0.094")	(0.094'')	(2.008'')	(2.619'')	3 bags per case



Part No. TX805	Head Width	Head Thickness	Head Length	Handle Width	Handle Thickness	Handle Length	Total Swab Length	Packaging
	19.7mm (0.776'')	9.4mm (0.370'')	26.4mm (1.039'')	6.6mm (0.260'')	3.2mm (0.126'')	105.8mm (4.165'')	132.2mm (5.205'')	80 swabs per bag 5 bags per case
							- in a	

Swab Solvent Compatibility									
Head Substrate Handle Material									
Solvent Foam Polyester Polypropylene Stat-Rite®									
Acetic Acid U U									
Acetone S U U D									
Acetonitrile U U D									
Ammonia, Anhydrous S U									
Aniline S U									
Benzene S U U D									
Chloroform S U									
hlorophenol D D									
Chromic Acid D D									
yclohexane S U									
limethylformamide D S									
Dimethylsulfoxide D S									
thyl Acetate S U									
thyl Alcohol U U									
thyl Chloride S U									
prmic Acid D U									
exane U U U U									
ydrochloric Acid 48% D D U U									
ydrogen Peroxide D U									
popropyl Alcohol U U U U									
lethanol S U U U									
lethyl Ethyl Ketone S U									
lethylene Chloride S U U D									
itric Acid 100% D D U D									
itrobenzene S U									
erchloroethylene U U									
hosphoric Acid 85% D U									
odium Hydroxide 36% S U									
ulfuric Acid 98% D D									
oluene S U U D									
richloroethylene S U									
Vater U U U									
Kylene S U D umpatibility is determined at room temperature Koy: D Domrados: S Swolls: U U									

Compatibility is determined at room temperature.

Sterile Available All Texwipe swabs available in sterile upon request. Key: D-Degrades; S-Swells, U-Unaffected



| Page 13







www.texwipe.com



Texwipe Americas+1 336 996 7046Texwipe Europe+31 88 1307 410Texwipe Asia+63 049 530 3390Texwipe China+65 6468 9433

www.texwipe.com

Rev.01-09/15